**ACTIVE SYSTEMS &
ADVANCED SYSTEMS**

TECHNICAL MANUAL

www.orthocormedical.com

TABLE OF CONTENTS

Technical Description	2
Package Components & Accessories	3
Indications for Use	3
Contraindications	4
Safety Information	4
Warnings & Cautions	5
Warranty & Efficacy	6-7
Frequently Asked Questions (FAQs)	8-9

Read instructions prior to use.

U.S. Federal law restricts the sale of this device by or on the order of a licensed health care practitioner.

Manufactured for
OrthoCor Medical - Blaine MN 55449

U.S. and foreign patents issued.
3.7 V DC internally powered.

INTRODUCTION

ABOUT THIS MANUAL

This manual describes the OrthoCor Active System or OrthoCor Advanced System and provides instructions for using it. Please read the entire manual carefully before use. This manual contains information, contraindications, and warnings, which must be observed to assure safe operation and to maintain the OrthoSource in a safe condition. This system **MUST** be used in accordance with the instructions and warnings.

TECHNICAL DESCRIPTION

The OrthoCor Active System or OrthoCor Advanced System are portable, compact, light-weight, non-invasive systems designed to deliver Pulsed Electromagnetic Field (PEMF) therapy for treating post-operative pain and edema in superficial soft tissue. The system consists of a multiple-use therapy device, referenced as the 'Active or Advanced System,' a multiple-use, rechargeable PEMF module, referenced as the 'OrthoSource,' and single-use therapy pods, referenced as 'OrthoPods' which are sold separately. The Active and Advanced Systems house both the OrthoSource and OrthoDocks which are appropriately positioned to deliver therapy to the treatment site.

Before using the either Active or Advanced Systems, please

read all of the indications, contraindications, warnings, and cautions in the following sections.

PACKAGE COMPONENTS & ACCESSORIES

Your OrthoCor Active or OrthoCor Advanced System contains:

OrthoSource: Module that contains electronics necessary to deliver Pulsed Electromagnetic Field (PEMF) technology

OrthoDock: Where OrthoPods are inserted on the Active or Advanced System

Active/Advanced System: Multiple-use device designed to house both the OrthoSource and OrthoPods and direct therapy to the appropriate areas of the body.

AC Adapter/Charger: Micro-USB compatible AC charger used to charge the OrthoSource. 3.7 V DC internally powered.

These therapy accessories are sold separately:

OrthoPods: Single-use therapy pods

INDICATIONS & CONTRAINDICATIONS

INDICATIONS FOR USE

The OrthoCor Active and Advanced Systems are indicated for adjunctive use in the palliative treatment of post-operative pain and edema in superficial soft tissue.

CONTRAINDICATIONS

- Do not use on patients with a cardiac pacemaker, cardioverter defibrillator, neurostimulator, infusion pump or any implantable electronic device.
- Do not use on patients who have ANY implanted metallic leads (wires) or any type of wire coil implant, or any implanted system that may contain a lead (wire).
- Do not use on patients who are or may be pregnant.
- Do not use on patients with an open wound at the area of application.
- Do not use on patients under the age of 18 or individuals with open bone growth plates.

SAFETY INFORMATION

SAFETY INFORMATION

- This device has been manufactured and tested in accordance to IEC 60601-1 and 60601-1-2.
- There are no serviceable parts inside the OrthoCor Active & Advanced Systems.
- Energy from PEMF therapy can be transferred through an implanted system and can cause tissue damage at the location of the implanted electrodes, resulting in severe injury or death.
- The OrthoCor Active System and OrthoCor Advanced System may interfere with appliances that are designed to detect and amplify weak radio signals. This includes television, radios, citizen band radios, phones, and other consumer devices in these categories. The interference may be heard as a repetitive tick in an FM radio and may cause tearing of a TV picture.

WARNINGS

WARNING: This device should be used under the continued supervision of a licensed health care practitioner.

WARNING: Use carefully. Remove immediately if discomfort occurs.

WARNING: Precaution should be taken when using this device on the elderly or incapacitated persons.

WARNING: Precaution should be taken when using this device on patients with sensitive skin. If the patient experiences skin irritation, discontinue the use of the device and consult your licensed health care practitioner before using the OrthoCor Active or OrthoCor Advanced System again.

WARNING: Do not over-tighten the Active or Advanced System as this may restrict blood flow and result in injury.

CAUTIONS

CAUTION: Use this device only in the prescribed manner and for the prescribed diagnosis.

CAUTION: KEEP AWAY FROM WATER. Use in dry environment.

CAUTION: The device is not indicated for treatment of deep tissue such as internal organs.

CAUTION: Do not use while the OrthoCor Active or OrthoCor Advanced System is charging.

CAUTION: Do not machine wash the OrthoCor Active or OrthoCor Advanced System.

CAUTION: Do not crush the OrthoSource or the Active/Advanced System as this may damage the electronics.

CAUTION: NEVER heat the OrthoCor Active or OrthoCor Advanced System or OrthoPods in the microwave as the system could catch on fire.

CAUTION: Removing the OrthoPods will stop the treatment.

CAUTION: The AC adapter/charger must be disconnected from the OrthoSource for treatment to occur.

CAUTION: Keep out of reach of children and pets.

CAUTION: For external human use only.

WARRANTY & EFFICACY

WARRANTY

OrthoCor Medical has a six-month manufacturer's warranty on any defects. OrthoCor Medical does not accept returns or offer refunds, except for cases that fall within the six-month warranty policy.

SUPPORT

If you have questions regarding the use of this device or our medical therapy, please email, call, write or visit our website.

OrthoCor Medical

8611 W 35W Service Drive NE
Suite 180
Blaine, MN 55449

888-583-6268

customersupport@orthocormedical.com

www.orthocormedical.com

CLINICAL EFFICACY

The Andrews Institute for Education & Research conducted a physician initiated, double-blinded, placebo controlled study of the use of the OrthoCor Active System for the treatment of pain for osteoarthritis.

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©2023 All rights reserved.